DIFFRAC.EVA 한글메뉴얼

By Bruker AXS Korea

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
A base line at the area bottom

A dashed line at the half maximum value (shows **FWHM** and the chord)

A cross representing the gravity center position

A vertical arrow representing the maximum position

Item	Description
Angle	
Left/Right Angle	Ends of the computation: angles. The input is rounded to the closest data point in X while the input is retained in Y
Obs. Max	Angle corresponding to the maximum intensity
FWHM	Full Width at Half Maximum value
Chord Mid.	Middle of the chord drawn between the mean values of the crossing points used to determine FWHM
I. Breadth	Net area (in cps×scan units) divided by the net height (in cps). It is the breadth of a rectangle having the same net height and the same surface as the peak. Given in scan unit
Gravity C.	Third estimate for the peak location. It is the center of gravity of the net peak. It is the mean of each X position in the interval weighted by the net intensity. It is also given in d (Å) if the scan is 2θ
Area (cps×deg)	
Raw area	Computed with the trapeze method and given in cps×scan units (cps×degrees for angular scans)
Net area	
Intensity (cps)	
Left/Right Int.	Left/right height given in cps
Gross Int.	Gross intensity
Net Height	Self-explanatory
Scherrer evaluation	
Crystallite Size	Crystallite size in Angstroms (Å)
Use FWHM	Select this option if you want to use the FWHM for the calculation of the crystallite size
Use I. Breadth	Select this option if you want to use the Intensity Breadth for the calculation of the crystallite size
K =	Scherrer constant for the calculation of the crystallite size with the Scherrer formula usually taken as 1 or 0.89. The default value is 1
Instr. Width =	Instrumental FWHM (FWHM for a material that exhibits no broadening beyond the instrument contribution) used for the calculation of the crystallite size with the Scherrer formula. The default value is 0

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
Filter	Description
Skip Non Ambient	Reject the patterns measured at a non-ambient temperature.
Quality marks	Reject patterns with a defined quality mark.
Status	Reject patterns with a defined status from the search.
Colors	Limit the possible color(s) of the compound searched for.
Sources	Select the sources to use for the search.
Subfiles	Limit the subfiles of the database to search in.
Element # in formula	Define a minimum and/or maximum number of elements in the compound formula.
Density	Define a minimum and/or maximum density for the compound formula.

Mark	Description
Low precision	Low precision pattern
Indexed	Good quality pattern with indexed lines
Blank	Pattern not meeting the Star, Indexed, or Low-Precision criteria
Star (*)	High quality pattern
Calculated	Pattern computed from single-crystal structural parameters
Rietveld	Pattern resulting from a Rietveld refinement
Hypothetical	Pattern calculated theoretically from an isostructural compound
Prototyping	Quality given to patterns which structural data was assigned by an editorial action of this Linus Pauling File (not recovered from the primary literature). This quality mark is specific to PDF-4

Status	Description
Primary	Pattern recommended by the PDF editorial board
Alternate	Pattern of reasonable quality which has been flagged as a duplicate of a primary pattern
Deleted	Pattern which has an improved replacement (a primary pattern). Its use is not recommended by the PDF editorial committee

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Search / Match (scan) m1.raw #1

Rebuild Chemical + - a/e Chemical Filter #1 Database + - a/e Database Filter #1

Database PDF2002: 148379 - After Filters: 132929

Chemical Filter Database Filter Candidate List Selected Candidates

Index #	Name	Formula	%	Source	ID
1	Calcite	Ca C O3	90	PDF2002	PDF 7
2	Calcite, syn	Ca C O3	69	PDF2002	PDF 0
3	Calcite	Ca C O3	69	PDF2002	PDF 4
4	Calcite	Ca C O3	34	PDF2002	PDF 2
5	Calcium Carbonate	Ca C O3	30	PDF2002	PDF 8
6	Carlinite, syn	Tl2 S	34	PDF2002	PDF 2
7	Magnesium calcite, syn	(Mg0.03 Ca0.97) (C O3)	15	PDF2002	PDF 8
8	Calcite	Ca C O3	77	PDF2002	PDF 0
9	Strontium Germanium Oxide	Sr Ge O3	44	PDF2002	PDF 4
10	Aragonite	Ca C O3 / Ca O · C O2	10	PDF2002	PDF 0
11	Aragonite, syn	Ca C O3	11	PDF2002	PDF 0
12	Calcite	Ca C O3 / Ca O · C O2	16	PDF2002	PDF 0

Group Duplicates Matched 38359 / 132929 Candid...

Search / Match

Whole Range Subrange Criterion: 2: Neutral

Auto Match

Search / Match Search by Name Filter Lists

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Column Description

Index # : Rank number of the given pattern during the search run

FOM (Figure Of Merit) : the higher this number, the better the match. Patterns are sorted by decreasing figures of merit (FOM)

Mtc : Number of reference pattern lines matching the unknown in the displayed range

nM : Number of reference pattern lines not matching the unknown in the displayed range

% : Actual Y-Scale percentage

Source : Name of the original database

ID : Pattern number in the database

Q : Quality mark

S : Pattern status

I/Icor : I/Icor is the ratio I/I_{cor} between the intensities of the strongest line of the compound of interest and the strongest line of corundum, both measured from a scan made of a 50-50 mixture, as stored in the PDF database

Mineral : Mineral compound or not

Name : Compound name as written in the database

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Selection : Patterns selected in the candidate list. Only the Delete operation can be applied to this selection

Doc. Patterns : All the patterns of the current document

Search by Name : All the patterns from the current Search by Name results

Search/Match : All the patterns from the Search/Match results

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Search / Match (scan) m1.raw #1

Rebuild Chemical a/e Chemical Filter #1 Database a/e Database Filter #1

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

Search / Match (scan) m1.raw #1

Rebuild Chemical + - ale Chemical Filter #1 Database + - ale Database Filter #1

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Filter	Value	Candidates
<input checked="" type="checkbox"/> Skip Non Ambient		138335
<input checked="" type="checkbox"/> Quality Marks		
<input type="checkbox"/> Status		
<input checked="" type="checkbox"/> Colors		
<input checked="" type="checkbox"/> Sources		
<input type="checkbox"/> Subfiles		
<input checked="" type="checkbox"/> Element # in For...		
<input checked="" type="checkbox"/> Density		

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Title	Value
Database	PDF2002
Pattern Number	PDF 72-1937
Formula	Ca C O3
Mineral Name	Calcite
Name	Calcium Carbonate
Quality	Calculated
Status	Primary
I/Icor	3.31
System	Rhombohedral
Space Group	R-3c (167)
a	4.994
c	17.081
Volume	368.93
Molecular Weight	100.09
Z	6
Density (calculated)	2.702
F(N)	999.9 (0.0001,30)
Subfiles	Inorganic, Mineral, Data source structure databases
Comments	
ICSD Collection Code	020179
Temperature Factor	ATF
Remark From ICSD/CSD	REM M PDF 5-586
Structure	
Publication	Dokl. Akad. Nauk SSSR
Detail	volume 245, page 1099 (1979)
Authors	Borodin, V.L., Lutin, V.I., Ilyukhin, V.V., Belov, N.V.
Primary Reference	
Publication	Calculated from ICSD using POWD-12++

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Data Tree

Data	Description	Views
Document		
Settings	1 Chemical Filter - 1 Database Filter	
2Theta	1 Scan	2
m2.raw #1	18P100 (Coupled TwoTheta/Theta)	
Pattern List #15	1 Pattern	
DIF (m2.raw)	18P100	


- Delete
- Select Parent
- Tool**
 - Search / Match (pattern)
 - Search by Number
 - d x by
 - Residue
 - Make Sticks
 - Make Peaks
 - Auto-scale
 - User Database**
- Create
- File

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The 'User Database Creation' dialog box contains the following fields and values:


- Short Name: BRUKER
- Long Name: Bruker Korea data
- Prefix: BK (example: ABC)
- Digit groups: 24 (example: 25 meaning: XX-XXXXX)
- Example: BK 45-6789

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The 'User Database DIF (m2.raw)' dialog box has the following configuration:

- Description: Bibliography, Experimental, Cell, Peaks
- Name: 18P100
- Formula: (empty)
- Keywords: (empty)
- Color: (empty)
- Subset: Inorganic (checked), Mineral (unchecked)
- User Database: BRUKER (dropdown), Bruker Korea data (text)
- Pattern ID: BK 00-0002 (dropdown)
- Buttons: Update ID, Delete ID, Append ID
- Checkbox: and Document

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